

No. 500
ECONOMY ADULT AIRWAY
MANAGEMENT TRAINER

Thank you for buying an American made product.

Recommendations:

The Economy Adult Airway Management Trainer Part Number 500 accepts up to 8.00 mm. E.T. Tube oral and nasal.

Set Up:

Your airway management trainer is ready for use. NO SET UP is required.

Lubrication:

Use water based lubricants only.

A tube of lubricant is included with the purchase of your airway management trainer. Place a coating of lubricant on the E.T. tube before inserting it into the trainer. Use of a petroeulement-based lubricant will void your Simulaids warranty.

Cleaning:

Simulaids airway management trainers can be cleaned with mild soapy water. Rinse thoroughly with clean water and let dry before storing.

Part Numbers:

No. 500 Economy Adult Airway Management Trainer

No. 4051 Lungs with Stomach 3 pack

No. 0363 Lubricant

Warranty:

Simulaids warrants this product to be free from any defect in materials and/or workmanship for a period of three years from the date of purchase, as evidenced by the date of invoice when the product was shipped to the end user. This warranty expressly does not cover abuse, accidental or purposeful damage, or any form of modification to the product. Simulaids reserves the right to either repair or replace affected parts or the entire unit, at their sole discretion, after investigating and reviewing the actual product and the damage. In most instances, a digital photo of the product in question showing the damage will help qualify a product for return to the factory. At no time will any product be accepted at the plant without proper return authorization issued by Simulaids.

Freight and Shipping charges are the sole responsibility of the end user. No product will be received with shipping charges due. Any product considered for warranty work must be identified by serial number and invoice number from the agency through whom the product was purchased. Without this information the product will not receive a return authorization number as required above.


Repair:

Simulaids airway trainers can only be repaired by Simulaids factory. Any other attempt to repair the trainer will void the warranty.

Return Goods Policy:

Should it be necessary to return the trainer to Simulaids for any reason contact the customer service department to obtain a RETURN AUTHORIZATION NUMBER.

Merchandise without a RETURN AUTHORIZATION cannot be accepted. When requesting your RETURN AUTHORIZATION please have your invoice number, date of purchase, or the product serial number. Do not pack trainer being returned in any type of packing materials that includes print like newspaper.



Manufacturers of Training Manikins,
Casualty Simulation Kits,
Medical Training Devices
PO Box 1289 - 16 Simulaids Dr
Saugerties New York 12477
(845) -679-2475

Toll Free: (800) -431-4310 Fax: (845) -679-8996
www.simulaids.com