

Model 26 Cot

*Simple, Easy to use, Independent Leg Cot
with Shortening Features*



MODEL 26 COT

The Ferno Model 26 is a simple, easy to use two-Level, roll-in Cot. It has only two positions - fully elevated, and ground level. In the elevated position, the cot can be shortened by elevating the backrest and folding down the hinged, head-end frame. Its legs fold up as it is rolled into the ambulance, and unfold as the cot is removed. An automatic locking system secures both the front and rear legs simultaneously as the leg sections unfold. Left-hand and right-hand control levers unlock the front and back legs respectively as the cot is prepared for loading.



FEATURES

- Simple to use roll-in cot with two levels of height adjustment: fully elevated or ground level
- Independent folding legs
- Fold-down drop frame shortens the length of the cot
- Lightweight tubular aluminum construction
- Swing-down sidearms
- Two stationary and two 360° swivel wheels
- Adjustable backrest angles from 0° to 75°
- Includes footrest, leg lock, and safety lock
- Available with bolster mattress, restraints, and shoulder harness for patient comfort and safety
- Now includes two roll-out wheels

SPECIFICATIONS

Length:	1920 mm
Width:	550 mm
Height:	Max: 800 mm
	Min: 610 mm
	Loading: 240 mm
Weight:	24.5 kg
Load Limit:	181 kg

ACCESSORIES & FASTENERS



360 Mattress Flat



460 Mattress Bolstered



430 Restraints



417-1 Restraints




513-10 IV Pole




514 Oxy-Clip




274 PacRac




175 plus 175-1




175 plus 175-2



175 plus 175-3



175 plus 175-4



175 plus 175-5