

# QUICK-RESPONSE COMPRESSED GAUZE

**SPECIFICATION** | 4.5" × 4.1yds stretched  
6-Ply, Cotton, White



## DESCRIPTION

RHINO RESCUE QUICK-RESPONSE COMPRESSED GAUZE provides a new idea of use, solves the problem that conventional compressed gauze cannot quickly identify and find out the position of the gauze during use, greatly increases the convenience of use and shortens the rescue time.

## INSTRUCTIONS FOR USE-A



Open from the top of the package, the gauze can be pulled out directly to fill the wound. If not all are used up, cut off what has been used and keep the rest in the package.

## INSTRUCTIONS FOR USE-B



Open from the bottom of the package, take out the entire compressed gauze, you can see the protruding part at one end of the gauze, you can quickly open the gauze by pulling out this protruding part.

